

VELP Scientifica Solutions for Stirring



MST
(Analog / Digital)

MICROSTIRRER

ESP

MULTISTIRRER 6
(Analog / Digital)

MULTISTIRRER 15
(Analog / Digital)

AMI

AMI4

Stirring speed rpm	up to 1500	up to 1100	up to 1100	from 80 to 1500	from 80 to 1500	up to 1100	up to 1100
Stirring volume L	up to 5	up to 5	up to 5	up to 2.4 (400 ml each)	up to 3.75 (250 ml each)	up to 5	up to 20
Timer	• (digital version only)			• (digital version only)	• (digital version only)		
Reverse Operation				• (digital version only)	• (digital version only)		
Dimensions (WxHxD) mm	120x50x145	120x48x128	160x33x230	230x51.5x370	230x51.5x370	150x55x270	600x55x270
	(in) (4.7x1.97x5.9)	(4.7x1.9x5.0)	(6.3x1.3x9.0)	(9.0x2.0x14.5)	(9.0x2.0x14.5)	(5.9x2.2x10.6)	(23.6x2.2x10.6)
Weight	Kg 0.6	0.55	0.9	1.75	2.1	1.2	4.0
	(lb) (1.3)	(1.1)	(2.0)	(3.8)	(4.6)	(2.6)	(8.8)
Power Supply	100÷240 V	100÷240 V	100÷240 V	100÷240 V	100÷240 V	100÷240 V	100÷240 V
Power	0.6 W	0.6 W	5 W	3.6 W	9 W	1.2 W	4.8 W



www.velp.com